

POLR3C Polyclonal Antibody

Description

Product type Primary Antibody

Code BT-AP07322

Host Rabbit

Isotype IgG

Size 20ul, 50ul, 100ul

Immunogen The antiserum was produced against synthesized peptide derived from human RPC3. AA range:191-240

Mol wt 60612

Species reactivity Human

Clonality Polyclonal

Recommended application WB, ELISA

Concentration 1 mg/ml

Full name POLR3C Antibody

Synonyms POLR3C; DNA-directed RNA polymerase III subunit RPC3; RNA polymerase III subunit C3; DNA-

directed RNA polymerase III subunit C; RNA polymerase III 62 kDa subunit; RPC62

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Specific core component of RNA polymerase III which synthesizes small RNAs, such as 5S rRNA and tRNAs. May direct with other members of the subcomplex RNA Pol III binding to the TFIIIB-DNA complex via the interactions between TFIIIB and POLR3F. May be involved either in the recruitment and stabilization of the subcomplex within RNA polymerase III, or in stimulating catalytic functions of other subunits during initiation. Plays a key role in sensing and limiting infection by intracellular bacteria and DNA viruses. Acts as nuclear and cytosolic DNA sensor involved in innate immune response. Can sense non-self dsDNA that serves as template for transcription into dsRNA. The non-self RNA polymerase III transcripts, such as Epstein-Barr virus-encoded RNAs (EBERs) induce type I interferon and NF- Kappa-B through the RIG-I pathway.

Recommended Dilution

WB: 1: 500 - 1: 2000 ELISA: 1: 40000

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using POLR3C Polyclonal Antibody cells nucleus.

Western blot analysis of RPC3 Antibody. The lane on the right is blocked with the RPC3 peptide.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com